




CAUTION!

- Before using the device please read the operating instructions!
- The instructions are part of the product and contain important information on getting started with and using the device.
- Always keep the instructions in a safe place for reference!
- They must also be passed on in the event the device is transferred to a third party.

Explanation of symbols

	CE mark		Green Dot Symbol
	GOST-R Certificate		Recycling Symbol
	Follow Operating Instructions		Manufacturing LOT
	Warning		Manufacturer
	Food Contact		EC Representative

Warning!

1. Keep this breast pump away and out of reach of children.
2. Clean breast pump with a soft, clean and nonabrasive cloth, do not use solvents or detergents.
3. If you suffer any pain or discomfort while using the device, please stop to operate and take a rest.
4. Always sterilize and clean all the parts prior to use in order to ensure the best health for you and your baby.

Parts




Assemble all the parts

1. Please use the breast pump in a total comfort.
2. Wash your hands and wipe your breasts especially around the nipple area with a clean warm towel. Massage your breasts and gently squeeze out a little breast milk to insure that your milk flow is fine.



- a) Close the breast shield to your breast tightly.
- b) Operate with handle in a comfortable rhythm.
- c) Adjust handle and Breast shield as you feel more comfortable.



d) When finish, take off the base and the bottle from the three way separator

e) Finally assemble the bottle in order to be ready to feeding your baby or store it in the refrigerator

⚠ Warning: It is important to check that all parts are ready compared to the component diagram.

Disassemble all the parts

1. Diagram as show below:



a) Take of the base from bottle

b) Twist off bottle from the three way separator


c) Take out the valve



d) Take out breast Shield

e) Rotate the handle, silicone and tank clockwise or counter clockwise by 100 degrees

f) Take out the handle, Silicone and tank from the three way separator.



- g) Lift up the handle and de-touch it from the three way separator (1)
Disconnect the handle from the silicone suction
- h) Separate the silicone from the tank

2. Cleaning and maintenance: After cleaning all the parts with a warm soapy water, disinfect the same parts with a sterilizer or boiling water for 10 minutes.

3. After the cleaning process is completed, please dry and reassemble all the parts as illustrated above.

Features

- Designed for all the mothers who wish to breast feeding.
- Safe, food grade material, non-toxic and BPA Free.
- Easy to operate: rotate the handle by 100 degrees to guarantee the best comfort.



Bredmed Limited
www.bredmed.com



Obebé S.r.l.
Viale Bianca Maria 23
20122 Milano, Italy.